EYELA

Stirrer Mazela

ZZ-1000 · ZZ-1000S ZZ-1100 · ZZ-1100S ZZ-1200 · ZZ-1200S ZZ-1300 · ZZ-1300S ZZ-1010 ZZ-1110 ZZ-1210 ZZ-1310

Instruction Manual

This instruction manual is designed to use the product efficiently and safely with keeping its best performance.

Be sure to read "Safety precautions" IMPORTANT before use.

Please keep this manual in a place easily accessible to every user.

IMPORTANT SAFETY INFORMATION

1. Warning signal word

Using a flammable or combustible liquid as a coolant for this product may result in a personal injury or an accidents when the product is improperly handled. Also, operating the unit under a relatively high room temperature may decrease its performance or cause a malfunction due to its functions or characteristics.

Most of such troubles, however, can be prevented if you are well informed about them beforehand.

To ensure the safety, this manual defines the information on such matters as requiring particular care in the safety as follows in terms of the importance and risk and attaches the alert mark and signal word.

It is recommended to follow the instruction to ensure the safe use of the product.

Alert mark Signal word	Definition			
	Danger of death or severe injury is expected when handled improperly.			
	Wrong handling is assumed to cause the possibility of the death or heavy injury of the user.			
Caution	Wrong handling is assumed to cause the risk of injury of the operator or physical damages.			

We have undertaken thorough verification concerning the possible occurrence of risk in the course of use of the product, but prediction of all and every kind of risk is extremely difficult. Namely, cautions contained in this manual are not necessarily all of possible risks.

However, if the product is operated according to the procedure described in this manual, the safe operation and work is ensured. Be sure to pay utmost care during handling of the product to prevent accident or failure of the product.


2. Warning Display on the Product

For particularly important warning instructions, the warning label is provided to the product main body.

The labelling position is shown below.

When using the product, be sure to pay due attention to the description of the warning.

* If damaged and illegible, be sure to change the warning label to the new one. Send the request for the new label to us.



This instruction manual explains installation, operation, troubleshooting, maintenance and inspection, and discarding procedures for the products below:

Always read this manual before use to ensure familiarization of the product.

Table of contents

1. For safe operation	1
2. Product overview	
2-1 Application	2
2-2 Operating procedures	2
2-3 Specification	2
2-4 Performance (Reference data)	3
2-5 Names of parts	4
3. Names and functions of the control assembly	
3-1 Control panel	5
3-2 Safety, Alarm and Warning functions	6
3-3 About the protection function	7
3-3-1 Thermal protector (Overheat protection circuit)	7
3-3-2 Current restriction circuit	7
4. Installation	
4 - 1 Installation environment	8
4 - 2 Installation conditions	8
4 - 3 Connecting the utilities	9

Mazela ZZ series


<ZZ-1000 series > · ZZ-1000 · ZZ-1000S · ZZ-1100 · ZZ-1100S · ZZ-1200 · ZZ-1200S · ZZ-1300 · ZZ-1300S <ZZ-1010 series > · ZZ-1010 · ZZ-1210 · ZZ-1110 · ZZ-1310

5. Operation

5-1 Pr	eparing for operation	10
5-1-1	Installing the main unit	- 10
5-1-2	How to open/close the chuck cover	- 11
5-1-3	Preparing the stirrer blade and the stirrer bar	12
5-1-4	How to install and remove the stirrer bar -	- 13
5-1-5	How to replace collets	- 16
5-1-6	Air purge	- 17
5-1-7	Adjusting the position of the main unit	- 18
5-1-8	Connecting the power plug	- 18
5-2 Op	perating methods	- 18
6. Causes	and solutions of troubles	- 19
7. Mainten	ance and Inspection	
7-1 Cle	eaning and care of the product	20
7-2 Ad	ljusting the chuck cover	- 20
8. Disposa	al of the product	- 21
9. After-sa	le Services	22
10. List of c	onsumable and replacement parts	- 23

Items contained in your packing

Be sure to confirm the type and quantity of parts before setting.



This product is not of an explosion-proof construction. Take extreme care for safety when operating this unit.

A Warning	Take extreme care for use of ignitable or flammable solutions flammable solutions will evaporate when stirred and heated and may ignite and explode if ignited at the motor assembly of the unit. When using such a solution, make a vent hole on the container and discharge vapor to the external through a tube or use the unit in a draft chamber and take extreme care.		
A Caution	Do not rotate the flexible stirrer in the left direction		
	Do not use a flexible stirrer specified for a drill chuck, which stirrer may cause the drill chuck to loose to allow the stirrer bar to fall and damage glassware when it is rotated left (anti-clockwise when the product is looked from above).		
▲ Caution	Close the chuck cover before using the unit		
	Never fail to close the chuck cover beforehand to prevent hair or clothing from being caught with the drill chuck or the collet chuck.		
▲ Caution	Check that the free clamp is secure before using		
	Operating with the loose free clamp will make the stirrer bar flatter and damage glassware. Be sure to securely tighten the knob bolt and the main unit is held in place before attempting operation.		
▲ Caution	Use the specified stirrer blades only		
	Never use a stirrer blade other than those specified by us (in particular, larger ones), which might cause a malfunction of the stirrer.		
	※Refer to the "10. List of consumable or replacement parts" on P.23 for stirrer blades recommended by us.		
▲ Caution	Do not repeat switching ON/OFF in a short time (within 15 seconds) Switching power ON/OFF in a short time (within 15 seconds) might cause a malfunction of the unit. Be sure to check the indication has completely disappeared before switching ON again to restart the unit.		

2-1 Application

WARNING

Never attempt to modify the product. Operate the unit for the specified purpose only.

An electric shock or a malfunction may result if the product is modified or used for any purposes other than that specified.

2-2 Operating procedures

You can set a speed according to the application. Turn the speed control knob on the control panel to set to a speed you want.

2-3 Specification

This product is a stirrer that supports mixing, solving or reacting of reagents.

You can change the speed during operation or stop operation when necessary. See "5-2 Controlling methods" on P.18 for details.

Model			ZZ-1000	ZZ-1000S	ZZ-1010
Туре			Standard		With reverse function
Pe	Speed range	rpm		50~3000	
erfor	May targua	N • m/rpm	0.098/3000		
erm	Speed range Max. torque	kgf • cm/rpm		1/3000	
	Motor output			DC brushless motor 35W	
Standard	Rated torque			0.098N • m(1kgf • cm)/3000rpr	n
dar	Shaft dia. [m	m]	φ1	~ 8	φ8(standard)
d.	Chuck assy.		Drill	chuck	Collet chuck
т	Rotation control		Speed feedback control		
Î	$\hat{\Xi}$ Speed setting and indication		Volume setting & digital indication: 4 digit		
Function	Safety function	on	Thermal protector (overhe	eat protection circuit), current res	triction circuit, chuck cover
	Included fund	sluded function Slow up/down, free clamp			
Oper	ating env. ten	np. range	5~35°C		
Ext. size[mm]		130W×155D×158H			
		1 * 2	$(130W \times 172D \times 173H)$		
Mass	5			Approx. 2.5kg	
Powe	er input & sou	rce voltage		AC100V 50/60Hz • 2A 200VA	

			•	1		
Model			ZZ-1100	ZZ-1100S	ZZ-1110	
Туре			Standard		With reverse function	
Pe	Speed range	rpm		20~1200		
ance	Mary tanana	N • m/rpm		0.245/1200		
e m	Speed range Max. torque	kgf • cm/rpm		2.5/1200		
	Motor output			DC brushless motor 35W		
Standard	Rated torque			0.098N • m(1kgf • cm)/3000rpm	n	
dar	Shaft dia. [m	m]	φ1	~ 10	φ8(standard)	
ď	Chuck assy.		Drill	chuck	Collet chuck	
ч	Rotation control Speed setting and indication Safety function			Speed feedback control		
l			Volume setting & digital indication: 4 digit		4 digit	
tio	Safety function	ety function Thermal protector (overheat protection circuit), current restriction circuit, chuck co		triction circuit, chuck cover		
B	Included fund	ction		Slow up/down, free clamp		
Oper	ating env. ten	np. range		5∼35 °C		
Ext. size[mm]		1 × 2	130W×155D×158H			
		1 * 2		$(130W \times 172D \times 173H)$		
Mass	5			Approx.2.5kg		
Powe	er input & sou	rce voltage		AC100V 50/60Hz • 2A 200VA	Δ	

X1 Sizes in () include protrusions.

X2 Depth when a free clamp is installed is 362[mm].

* The performance shown is at a room temperature of 20°C, rated source voltage, 50Hz and no-load.

Model			ZZ-1200	ZZ-1200S	ZZ-1210
Туре			Stand	lard	With reverse function
ге	Speed range	rpm		10~600	
ance	Speed range Max. torque	N • m/rpm		0.49/600	
e m	wax. torque	kgf • cm/rpm		5/600	
S	Motor output			DC brushless motor 35W	
tan	Rated torque		0	0.098N • m(1kgf • cm)/3000rpi	n
tandard	Shaft dia. [mi	n]	φ1~	-10	φ8(standard
d	Chuck assy.		Drill c	huck	Collet chuck
Ч	Rotation control		Speed feedback control		
Function	E Speed setting and indication		Volume setting & digital indication: 4 digit		
tio	Safety function	on	Thermal protector (overhea	t protection circuit), current res	triction circuit, chuck cover
'n	Included func	luded function Slow up/down, free clamp			
Ope	rating env. ten	np. range		5∼35°C	
Evt	size[mm] 💥	1 ** 2		$130W \times 155D \times 158H$	
Ext. size[mm] *1 *2		×1 × Z	(130W×172D×173H)		
Mass				Approx. 2.5kg	
Power input & source voltage		irce voltage	_	AC100V 50/60Hz • 2A 200VA	A

Model			ZZ-1300	ZZ-1300S	ZZ-1310
Туре			Standard		With reverse function
	Speed range	rpm		5~300	
Perform		N • m/rpm		0.98/300	
em	Max. torque	kgf • cm/rpm		10/300	
s	Motor output	t		DC brushless motor 35W	
tan	Rated torque			0.098N • m(1kgf • cm)/3000rpn	n
Standard	Shaft dia. [m	m]	φ1	~ 10	φ8(standard
d	Chuck assy.		Drill	chuck	Collet chuck
Ч	Rotation con	trol		Speed feedback control	
Function	Speed setting and indication		Volume setting & digital indication: 4 digit		
tio	Safety functi	on	Thermal protector (overhe	eat protection circuit), current rest	triction circuit, chuck cover
n	Included fun	ction		Slow up/down, free clamp	
Opera	ting env. tem	p. range		5∼35℃	
Evt a	ize[mm] 💥	1 * 2	130W×155D×158H		
Ext. size[mm]		1 * 2	(130W×172D×173H)		
Mass				Approx. 2.5kg	
Powe	r input & sou	rce voltage		AC100V 50/60Hz • 2A 200VA	

%1 Sizes in () include protrusions.
%2 Depth when a free clamp is installed is 362[mm].
% The performance shown is at a room temperature of 20°C, rated source voltage, 50Hz and no-load.


<accessories></accessories>	
ZZ-1000,1100,1200,1300	Clamp x 1, chuck handle x 1
ZZ-1000S,1100S,1200S,1300S	Clamp x 1, chuck handle x 1 Stirrer bar x 1 (ϕ 8 × 495mm, SUS), stirrer blade of marine type x 1 (ϕ 50mm, SUS)
ZZ-1010,1110,1210,1310	Clamp x 1
2-4 Performance	
(Reference 3000 data) 2000 [rpm] 1000-	Speed — Torque characteristics — Model ZZ-1000,1000S,1010 — Model ZZ-1100,1100S,1110 — Model ZZ-1200,1200S,1210 — Model ZZ-1300,1300S,1310
_	
L	
0	0.1 0.2 0.3 0.4 0.5 0.6 0.7 0.8 0.9 1.0
	Torque[N · m]

2-5 Names of parts

ZZ-1000	ZZ-1000S	ZZ-1100	ZZ-1100S
ZZ-1200	ZZ-1200S	ZZ-1300	ZZ-1300S

< Front >

< Back >




X1 See "3-1 Control panel" on P.5 for details of the control panel.

3 Names and functions of the control assembly


3-1 Control panel

ZZ-1000	ZZ-1000S	ZZ-1100	ZZ-1100S
ZZ-1200	ZZ-1200S	ZZ-1300	ZZ-1300S



No.	Name	Description of the functions
1	Digital display	Indicates the current speed. Indicates an alarm No. when a malfunction occurred to the unit.
2	Speed control knob	You can change the speed by turning this knob. The further you turn this clockwise, the higher the speed of the unit will be.

ZZ-1010 ZZ-1110 ZZ-1120 ZZ-1130



No.	Name	Description of the functions	
1	Digital display	Indicates the current speed. Indicates an alarm No. when a malfunction occurred to the unit.	
2	[Fow/Rev] key Toggles the rotation direction [Left $\leftarrow \rightarrow$ Right].		
3	Speed control knob	You can change the speed by turning this knob. The further you turn this clockwise, the higher the speed of the unit will be.	

3-2 Safety, Alarm and Warning functions

This product has the following safety functions.

When any abnormality occurs, take appropriate measures by referring to the section "Causes and solutions of troubles" on P.19.

Safety function

Safety device	Operation	Causes	Solutions
Fuse	Melts and turns power off.		Contact your dealer or the nearest service center. %The device is installed in the board and the user is not allowed to replace it.
Current restriction circuit (overload)	Rotation indication flashes. (Operation continues at a lower speed)	The motor was subjected to a load beyond the rating.	Reduce load.
Lock protection	Stops the motor.	The motor output axis was locked (seized) from overloading.	Reduce load or eliminate the cause of lock.
Thermal protector (Overheat protection circuit)	Shuts off current to the motor.	The motor overheated during continued operation under overload.	Allow the unit to cool down (about 30 minutes).

Alarm function

** Alarm indication[A-**] : Appears to notify the operator of a tentative malfunction.

	-	
	Operation	Stops rotation (motor).
	operation	Digital displays indicates an alarm.
	Causes	Continued operation under overload resulted in overheat of the motor or
	Causes	the control board beyond the specified temperature.
[A-18] Overheat alarm	Solutions	After the unit has cooled down (about 30 minutes), restart it (turn power
	Solutions	on again).
		Stars actation (mater)
	Operation	Stops rotation (motor).
/-/ - _/ '-/	- r	Digital displays indicates an alarm.
	0	Current restriction circuit(overload) caused the motor speed down below
	Causes	10rpm for two or more seconds consecutively.
[A-19] Overcurrent alarm	Solutions	After reducing load, restart it (turn power on again).
	Operation	Stops rotation (motor).
	operation	Digital displays indicates an alarm.
/-/ - _' /_/	Cause	Source voltage to the unit is high.
		Turn the power switch OFF and check that the source voltage is within
[A-20] Overvoltage alarm	Solution	the commercial source voltage (or AC100V) $\pm 10\%$.
[]		
		Stops rotation (motor).
	Operation	Digital displays indicates an alarm.
└─┘└─┘└━┘└─┘	Cause	Source voltage to the unit is low.
	Solution	Turn the power switch OFF and check that the source voltage is within
[A-21] Low voltage alarm	Solution	the commercial source voltage (or AC100V) $\pm 10\%$.

Alarm function

※ Alarm indication[F-**] : Appears when an abnormality that is detected from a malfunction of a part requiring repair.

Operation	Stops rotation (motor). Digital displays indicates an alarm.			
Cause	Abnormality (disconnection, malfunction) occurred to the rotation detection sensor.			
Solution	Replacement of the motor (rotation motor) is necessary.			

[F-4] Rotation sensor error alarm

3-3 About the protection function

3-3-1 Thermal protector (Overheat protection circuit)

- (1) Continued overload operation raises the motor or the control board temperature and when the temperature detection element on the board detects that the temperature has reached the threshold, the protection circuit will activate and shuts off current to the motor. The display will indicate [A-18](Overheat alarm).
- (2) To recover, turn the power switch OFF once and wait until the motor or the control board temperature cools down.
- (3) Operation may be resumed when the motor temperature drops below the specification. (Rough standard: About 30 minutes) When you operate the unit again, lower the load than the setting when the overheat protection function was activated.

3-3-2 Current restriction circuit

- (1) When load exceeded the permissible range from its fluctuation during stirring, the circuit automatically lower the stirring speed to prevent overloading. During deceleration from overload, the current speed flashes in the display.
- (2) While the speed is flashing, you can operate the speed control knob to decelerate the speed and stop flashing.
- (3) Current restriction circuit keeps the stirring speed at 10rpm or below and overload does not improve, the display will show [A-19](Overcurrent alarm) and the motor will stop.
- (4) To cancel the alarm, turn the power switch OFF once, reduce load and then restart the unit.

Precise Torque (load) detected is not necessary precise and overload might be detected and the speed might flash even during normal operation.

4 Installation

4-1 Installation environment

WARNING

Do not install the product in a potentially hazardous location.

This products is not designed with explosionproof structure. Using the unit in a potentially hazardous location may cause a fire.


Select the installation site that meets the following conditions.

- Securely install using a special stand base (model ZS) or the existing angle on a level and stable surface.
- No inflammable solid or liquid or gas around the unit.
- Ambient temperature must be from 5 to $35 \,^{\circ}\text{C}$.
- No dew condensation
- Lesser humidity and no dripping on the unit.
- Lesser dust
- No direct sun light
- Well-ventilated

CAUTION

Install the product in the appropriate installation site.

Installing the product in he inappropriate location may accelerate the deterioration of the product or worsen the functions.




4-2 Installation conditions



Assure certain space around the unit.

Be sure to assure space between the unit and the walls and the ceiling as shown in the figure to maintain the performance of the product.



4-3. Connecting the utilities

WARNING

Confirm the voltage, phase, capacity, and the type of receptacle of power supply.

Wrong connection of power supply may cause fire or electric shock

▲ WARNING

Be sure to connect the earth wire.

Failure to connect the earth wire may cause an electrical shock.

- (1) Check the voltage, phase and capacity of the power source. Required power source is as shown in right table.
- (2) Check the receptacle of installation place.

Prepare a receptacle g with earth

- Do not connect the power plug in this stage.
- Check that the sheath of the power cord is not damaged.A damaged cord may cause an electrical shock.
- * Do not use a branching socket or a power strip for connection to the power supply.

Power cord specifications (common for all models)

	Model	Ca	able	Section area
		Length	Thickness (O.D.)	of cable(AWG)
	ZZ-1000 ZZ-1010 etc.	Approx.2.0m	Approx.7.4mm	0.75mm ² (AWG18)

A WARNING

Do not use the branching socket or table tap.


It may cause burnout of burn out of cables or a fire from overcurrent.

WARNING

Connect the earth wire correctly.

Never connect the earth wire to the gas pipe or water pipe to prevent electric shock.

Power supply for connection		
Voltage Current		
AC100V	15A	



5 Operation

5-1 Preparing for operation

5-1-1 Installing the main unit

- (1) Loosen the knob bolt on the stand base, insert the pole and retighten the bolt to securely fix the pole.
 - * The stand base and the pole are not included with the product.

<Stand base set>

Code No	Model	Pole size [mm]
158710	ZS	$\varphi 22 imes 780$
203190	ZS-2	φ22×1200

A Caution

Take sufficient care during transportation

The stand base is heavy.

Hold both the pole and the base for transportation If you hold only the pole, the base may drop and cause a personal injury.

Marning

Securely fix the clamp


Loose tightening of the clamp may allow the stirrer to vibrate and fall off.


(2) Install the clamp included with the product to the pole of the stand base assembled in (1) or to the existing angle and the securely tighten the knob bolt of the clamp to fix.

Possible angle diameter of the clamp : O.D.8 \sim 22mm

(3) Insert the arbor of the main unit through the clamp and then securely tighten with the knob bolt to fix.

Support the main unit with a hand during fixing work to prevent it from falling.





A Caution

Do not open/close the chuck cover during operation

If the chuck cover opens or closes while the chuck is rotating, touching it might cause a personal injury.


Be sure to check that rotation has stopped before opening the cover for replacing the shaft or other works.

- * Be sure to open or close the chuck cover only after turning power of the unit OFF.
- Using fingers, open the chuck cover slightly.
 * The open-preventive claw at the front of the cover is disengaged and the lock is released.
- (2) Open the right and left covers widely until the protrusion on the back enters into the center groove in the cover.
 - * When it is completely in the groove, the cover will be kept open enabling easier work.

(3) Push the whole chuck cover assembly backward.

- (4) Follow the procedures (1) to (3) above in the reversed order to put the cover over the chuck.
 - * You will hear a click when the chuck cover claw is engaged.

If no click is heard, it is not engaged and the chuck cover may open during operation. In this case, open it once and then close the chuck cover again.



Push up backward

5-1-3 Preparing the stirrer blade and the stirrer bar

※ "P.23 10.List of consumable and replacement parts" shows a list of stirrer blades and stirrer bars recommended by us.

Marning

Use the specified stirrer blades only


Never use a stirrer blade other than those specified by us (in particular, larger ones), which might cause a malfunction of the stirrer.

< Stirrer bar for stirrer blade >

(1) Turn the fixing screw on the side of the stirrer blade anticlockwise (left) to loosen.

XYou do not have to remove the screw.

- %The fixing screw is a very small part. Take care not to lose it by loosening it overly.
- (2) Insert the stirrer blade through the stirrer bar.
- (3)Loosen the fixing screw on the side of the stirrer blade clockwise (right) to tighten to fix it to the stirrer bar.



< Stirrer bar with screw>

% This is an accessory to the model with a stirrer bar and a stirrer blade 【ZZ-1 \square 00S】.

A Caution

Do not use a stirrer bar with a screw on models with the reverse function

The stirrer blade may drop off on reversal causing damages to containers or a personal injury.

- (1)Insert the screw through the hole of the stirrer blade.
- (2)Turn the screw clockwise (right) into the screw hole at the tip of the stirrer bar and sand the stirrer blade to fix.
- * Align the hole shape of the stirrer blade with the screw shape before fixing.

5-1-4 How to install and remove the stirrer bar

ZZ-1000	ZZ-1000S	ZZ-1100	ZZ-1100S
ZZ-1200	ZZ-1200S	ZZ-1300	ZZ-1300S

A Caution

Securely tighten the drill chuck

Securely install the stirrer bar to the chuck and fix with the chuck handle.

If it is not secured, the stirrer bar may slip or drop off and damage the reagent containers or cause a personal injury.

< How to install (drill chuck)>

- (1) Turn the chuck clockwise (right) stirrer bar with a hand while holding it with a hand to prevent it from coming off-center and lightly fix the stirrer bar.
- ※Insert the stirrer bar at least to the middle point of the whole length. If it is unfeasible, insert the bar so that it is protruding above the top of the unit.
- (2) You can install without off-centering by turning the stirrer bar right and left for 2 to 3 times when the chuck is lightly tightened.
- (3) Insert the chuck handle into the hole with a "〇" mark and securely tighten.


< How to remove (drill chuck)>


- Insert the chuck handle into the hole with a "○" mark while holding it with a hand and turn the chuck anticlockwise (left).
- (2) When the chuck is loosened and the claw has opened, gently pull out the stirrer bar downward.

A Caution

Do not use a warp stirrer bar or install it off-centered.

Operating the unit with a warp stirrer bar or installed off-centered may cause the unit vibrate, damage the reagent containers or the unit to fall off.





A Caution

Take care for the O.D. of the stirrer bar to use

Collet chuck only supports stirrer bars with a specified O.D. and cannot tighten other bars. Be sure to check the O.D. of the stirrer bar before use. (Standard O.D. : φ 8)

$\underline{\wedge} \, \textbf{Caution}$

Securely tighten the collet chuck


Use a stirrer bar of the same as the collet and turn it to fix securely.

If it is not fixed securely, it may slip, drop off, damage the reagent containers and cause a personal injury.

$\underline{\wedge}$ Caution

Do not use a stirrer bar with a screw on models with the reverse function

Use of a stirrer bar with a screw shown below may cause the stirrer blade to drop off on reversal causing damages to containers or a personal injury with broken pieces.



A Caution

Check the collet size

When the size of the stirrer bar and of the collet, the bar may not be fixed securely and it may slip, drop off, damage containers and their pieces may cause a personal injury.

A Caution

Do not use a warp stirrer bar or install it off-centered.

Operating the unit with a warp stirrer bar or installed off-centered may cause the unit vibrate, damage the reagent containers or the unit to fall off.

A Caution

Check for flaws or dirt on the collet or the stirrer bar

Use of a collet or a stirrer bar with an apparent flaw or dirt on its surface may prevent the collet from tightening evenly and lead to slip of the stirrer bar. Continuing operation with that state may lead to damages to parts or devices, contamination of metal chips into reagents, or off-centering of the stirrer bar.


< How to install (Collet chuck) >

 Check for flaws or dirt on the surface of the stirrer bar.
 Dirt or flaws may prevent it from properly installing to the collet chuck.

- (1) Loosen the collet chuck by turning the nut cover anti-clockwise (left) while holding the holder cover.
- (2) After inserting the stirrer bar slowly, loosen the nut cover by turning it clockwise (right) to fix the stirrer bar.
 - X You may not able to insert the stirrer bar even if you loosen the collet chuck.

This is because the collet inside is engaged with the holder cover, which can be released by slowly inserting the stirrer bar from above until it lightly hits and then lightly hit the collet to disengage it.

Insert the stirrer bar at least to the middle point of the whole length. If it is unfeasible, insert the bar so that it is protruding above the top of the unit.




< How to remove (Collet chuck) >

- (1) Loosen the collet chuck by turning the nut cover anticlockwise (left) while holding the holder cover.
- (2) When the chuck is loosened, gently pull off the stirrer bar downward.
 - You may not able to pull out the stirrer bar even if you loosen the nut cover collet chuck. This is because the collet inside is engaged with the holder cover.

In this case, loosen the nut cover further and repeat pulling the stirrer bar several times to disengage the collet.

If this does not work, remove the collet by referring to the section "P.16 5-1-5 How to replace the collet".



5-1-5 How to replace collets

ZZ-1010 ZZ-1110 ZZ-1210 ZZ-1310

A Caution

Check for flaws or dirt on the collet

Do not use a collet or a stirrer bar with an apparent flaw or dirt on its surface which may prevent the collet from tightening evenly and lead to slip of the stirrer bar.

Continuing operation with that state may lead to damages to parts or devices, contamination of metal chips into reagents, or off-centering of the stirrer bar.

- (1) Loosen the nut cover by turning anti-clockwise until the drop-off preventive screw on the side is exposed while holding the holder cover.
- (2) Loosen the drop-off preventive screw with a precision screwdriver.
 - XYou do not have to remove the screw.

*Drop-off preventive screw is a very small part. Take care not to lose it by loosening it excessively.

- %Precision screwdriver is not included. The user shall prepare them.
- (3) Loosen the nut cover further and gently remove the cover and the collet.
- (4) Take out the collet out of the nut cover and replace with a correct collet.


Code No	Product name	Supported shaft dia.(mm)
256100	Colletø4	φ3~4
256110	Colletø6	φ5~6
256120	Colletø8(standard)	φ7~8
256130	Colletq10	φ9~10

<List of replacement collet>

- ※ Collet I.D. corresponds to product names. Check the size before replacement or check if the size is correct by inserting it through the stirrer bar.
- (5) After inserting the collet into the nut cover, turn the cover clockwise (right) until just before the drop-off preventive screw is hidden to fix it to the holder cover.
- (6) Securely tighten the drop-off preventive screw with a precision screwdriver and turn the nut cover clockwise (right) only slightly to tighten.
 - Excessively tighten the nut cover without a stirrer bar will make the collect engaged with the holder cover, which might prevent the collet from disengaging and you cannot insert a stirrer bar by loosening the nut cover.

In this case, remove the nut cover out of the holder cover and pull out the collet.

Or by slowly inserting the stirrer bar from above until it lightly hits and then lightly hit the collet to disengage it.



A Caution

Avoid applying stress on the tube

Applying stresses to the tube such as tension, twist or a sharp bend may cause damages or coming-off of it.

This unit has an air purge port on the back. You can prevent corrosive gas from entering into the unit by introducing pressurized air at 0.01MPa or below at the air purge port. Use this function as necessary.

<Tube to use>

Use a pneumatic piping tube having an O.D. of 6mm. Also check that the tube meets the following conditions.

- The cut surface is square.
- There is not flaw on the outer surface.
- The section is not of an oval shape.

Also check precautions for the tube before use.

<Installing a tube>

- (1)When a plug is attached to the air purge port, pull out the plug while pressing the open ring.
- X Store the removed plug securely since it will be used when air purge is not performed.
- (2) Insert the tube into the air purge port until it stops.
- (3)After attachment, lightly pull the tube to make sure it will not come off.
- (4)Introduce pressurized air of 0.01MPa or less. Use high quality air free of oils or humidity.

<Removing a tube>

- (1) Stop air supply.
- (2) Pull out the tube while pressing the open ring of the air purge port.
- *Leave the air purge port inserted when air purge is not carried out.

A Caution

Keep the introduction pressure to 0.01MPa (0.1kgf/cm²) or lower of a gauge pressure


Introducing highly pressurized air might damage the stirrer.

A Caution

Use high quality air

Use of air that contains oil, humidity or corrosive gas may cause a malfunction.

Use a filter or a dryer to remove any foreign objects in air to introduce.



5-1-7 Adjusting the position of the main unit

Loosening the arbor knob bolt will also loosen the main unit and you can freely adjust its angle and the position. After adjustment, tighten the bolt to securely fix the main unit.

5-1-8 Connecting the power plug

Make sure that the power switch is OFF, turn the Speed control knob to the left most and connect the power plug to an outlet.

5-2 Operating methods

A Caution

Stop operation when an abnormality occurs

When an abnormality occurs, immediately turn the power switch OFF and refer to the section "P.19 Causes and solutions of troubles".

- Make sure that the chuck cover is closed and the speed control knob is at the leftmost end before starting operation.
- (1) Turn the power switch ON. The display will come on.
 - Note that if the speed control knob is not at the leftmost position, turning the power switch ON will start rotation of the stirrer blade.
- (2) While observing the display, slowly turn the speed control knob until the speed you want is obtained. The more you turn it clockwise, the faster the rotation will be.

%For a model with the reverse function $[ZZ-1\Box 10]$


Each time you press the [Fow/Rev], the rotation direction will change. (Left \rightarrow Right \rightarrow Left \rightarrow ...)

Set a rotation direction you want before starting operation.

- Be sure to reduce the speed before changing the direction.
 Suddenly changing the direction during stirring (rotation), load will be applied to the gears and cause wear.
- (3) TO finish, turn the speed control knob leftmost and stop operation then turn the power switch OFF.
 - Make sure that the display has turned off completely before turning power on again.

Procedure after operation

If you are not going to use the unit for an extended period of time, turn the switch OFF and remove the power plug off the outlet.



A Caution

Fix the stirrer container before use

Fix the container since high viscosity material may be used or some material will increase their viscosity during stirring causing the container to move.

A Caution

Keep the chuck cover closed during operation

Never fail to close the chuck cover beforehand to prevent hair, clothing or hands from being caught with the drill chuck or the collet chuck.

A Caution

Do not move or transport the unit during operation

Be sure to stop operation before moving or transporting the unit when touching the rotating assembly might cause a personal injury.

A Caution

Do not keep operating at the maximum performance

Stirring under a higher load at a high temperature might shorten the use life of the stirrer. Assure sufficient performance margin, in particular, at a higher temperature.

6 Causes and solutions of troubles Contact your dealer or the nearest service center for troubles not listed here.

Symptom	Causes	Solution	
	The power plug has come off from the outlet. Or it is not securely inserted.	Turn the power switch OFF and insert it to an outlet.	
Turning the neuror quitch ON will	Power is not supplied.	Turn the ELB on the distribution board ON.	
Turning the power switch ON will not show anything on the display.	The fuse is blown. (It is installed inside the product and cannot be checked from outside.)	Stop operation immediately and contact your dealer of the nearest service center.	
	The board is defective.	service center.	
Does not rotate even if the speed control knob is turned. (Display shows [A-18])	The thermal protector (Overheat protection circuit) activated from overload.	Turn power OFF once and allow the motor to cool down (about 30 minutes) First reduce load and then restart the unit.	
Does not rotate even if the speed control knob is turned(Display shows [A-19])	The motor spindle locked from overload.	Turn the power OFF once, reduce load and then restart the unit.	
Does not rotate even if the speed	Gear for conveying rotation from the motor is damaged.		
control knob is turned(rotation lamp on)	The speed control volume is malfunctioning.	Stop operation immediately and	
(Display shows the speed)	The board is defective.	contact your dealer of the nearest service center.	
Does not rotate even if the speed control knob is turned(Display shows [F-4])	The motor rotation sensor is malfunctioning.		
Speed dropped suddenly.	Current restriction circuit(overload) activated.	Turn the power OFF once, reduce load and then restart the unit.	
Abnormal sound is heard from the unit.	Gear for conveying rotation from the motor is damaged.	Stop operation immediately and contact your dealer of the nearest	
The unit is vibrating excessively.	The motor is malfunctioning.	service center.	
Display shows [A-20]	Source voltage to the unit is high.	Make sure that the source voltage is within the commercial power	
Display shows [A-21]	Source voltage to the unit is low.	supply(or AC100V) $\pm 10\%$.	
The stirrer bar cannot be removed by loosening the collet chuck.	The collet is engaged and cannot be loosened.	See the section "How to remove (Collet chuck)" on P.15.	

7-1 Cleaning and care of the product

A Warning

Never attempt to disassembly the product.

The unit contains parts with high voltage applied or may become hot, and disassembly may cause an electrical shock or an injury.

A Caution

Do not leave the product in an undesirable atmosphere

Leaving the product at a place where acid or organic solvent density is high will cause discoloration of the finish of the product or peeling-off of the paint. It also accelerate corrosion of the board resulting in a malfunction.

- (1) Before any maintenance work, turn the switch OFF and remove the power plug from the outlet.
- (2) Use a moistened and well wriggled soft cloth for cleaning. Use milt detergent for stubborn dirt and completely wipe remaining detergent after cleaning.

7-2 Adjusting the chuck cover

Repeated open/close of the chuck cover will make the connection to the main unit loose and may prevent fixing at a position you want.

If you notice loosening when you lift the chuck cover for replacing the stirrer bar, or the cover falls below the original position, turn the connection screw or the nut clockwise (right) to tighten and adjust the cover.

Marning

Do not clean and care the product while the unit is still hot.

Touching hot surface may burn your hands.

A Caution


Use a correct method and items for cleaning or caring the product.

When cleaning or maintaining the product, never splash water to the exterior or the inside directly, do not put any foreign materials and never use the cleanser, thinner, oil, kerosene, acid, and equivalent. Otherwise, the user may suffer electric shock or damage to the product.

A Caution

Frequently remove dusts from the outlet

Operating the unit for a long time with connected to an outlet will cause dust to accumulate between the outlet and the power plug and may lead to a malfunction or a fire. (Tracking)



8 Disposal of the product

When disposing the product, please follow the instructions as below.

Component	Model	Weight	External dimensions(including protrusions)	Method for disposing
	ZZ-1000 ZZ-1000S ZZ-1010	Approx. 2.5kg	130W×155D×158H (130W×172D×173H)	
	ZZ-1100 ZZ-1100S ZZ-1110	Approx. 2.5kg	130W×155D×158H (130W×172D×173H)	Request the disposal operator for
Main unit	ZZ-1200 ZZ-1200S ZZ-1210	Approx. 2.5kg	$\begin{array}{c} 130W \times 155D \times 158H \\ (130W \times 172D \times 173H) \end{array} disposal. \end{array}$	
	ZZ-1300 ZZ-1300S ZZ-1310	Approx. 2.5kg	130W×155D×158H (130W×172D×173H)	

Main com	onents and	l disposal	instructions
main com		usposa	1130 000013

Component	Major	Major materials
	components	
	Housing	Die-cast aluminum
Main	Chuck cover	PC
Main	Free clamp	Die-cast aluminum, stainless steel, brass
	Driving assy.	Stainless steel, aluminum, POM, MC-901, S45C, ABS, cupper

XWhen disposing of the product, separate materials according to the table above before disposal.

9 After-sale Services

- 1. In case the product does not function satisfactorily, check first by referring to the page on troubleshooting to see if this is actually a trouble.
- 2. If the product remains unsatisfactory even after checking, contact the shop from which the user has purchased the product or the service center described in the manual and request repair.
- 3. Repair during the guarantee period will be made according to the guarantee stipulations.
- 4. After expiration of the guarantee period, the charged repair will be made at the customer's request.

Warranty Provisions

- 1 The warranty period of your product is 12 months from the date of purchase.
- 2 We will repair or replace your product free of charge for troubles occurred in your product during proper usage in the warranty period.
- 3 This warranty is limited to the product itself and shall not be construed to cover all and any consequential losses (operating losses, various expenses, etc.) from malfunctions or troubles of the product.
- 4 The provisions of this warranty are valid only in Japan. When you are going to indirectly export the product you purchased to a foreign country, the provisions shall be exempted from application provided that a certification be issued that the product is not subject to the export control regulations and, in such case, all liabilities in relation to the product shall be borne by the exporting party.
- 5 Repair shall be paid by the user even during the warranty period for the following cases.
 - a) The user did not return the customer card or made user registration at our HP within one month from the date of purchase.
 - b) This warranty card was not presented or user registration could not be confirmed at the time of request for repair.
 - c) When the dealer name is not sealed on the warranty card or the date of the purchase is not indicated on it.
 - d) Malfunctions or damages caused by handling that is against cautions in the manual or labels on the product, or transportation from the installation site after purchase, dropping during use, or shocks.
 - e) Malfunctions or damages caused by mishandling, unauthorized modification or repair by the user.
 - f) Malfunctions or damages caused by a fire, an earthquake, wind and flood, salt damage, lightening, or other Acts of God or external factor including the power supply.
 - g) Degrading of performance or malfunctions resulted from consumption of consumable parts or replacement of consumable parts.
- 6 Warranty conditions may differ depending on specific products. Check the warranty conditions stipulated in the warranty space in the operation manual.
- 7 Ask your dealer or the nearest sales office for repair for malfunctions after the warranty period.
- (In principle, the retention period of repair parts is five years after the end of production.)
- 8 Product warranty for products sold overseas by our overseas sales department shall be separately specified.

10 List of consumable and replacement parts

Stirrer blade

	No.	No. 1	No. 2	No. 3	No. 4	No. 5	No. 6	No. 7	No. 8	
	uter arance	Į.	E.	Ð	×	K	X	×	Z	
Name		Flat blade turbine	Curved flat blade turbine	Disc turbine	Fan turbine	Curved fan turbine	Angled fan turbine	Modified fan turbine	Marine blade	
Axis dia. [mm]	Blade dia. [m]	Code No.	Code No.	Code No.	Code No.	Code No.	Code No.	Code No.	Code No.	
	30	121820	122000	122180	122360	122540	122720	122900	123080	
	40	121830	122010	122190	122370	122550	122730	122910	123090	
	50	121840	122020	122200	122380	122560	122740	122920	123100	
	60	121850	122030	122210	122390	122570	122750	122930	123110	
8	70	121860	122040	122220	122400	122580	122760	122940	123120	
0	75	121870	122050	122230	122410	122590	122770	122950	123130	
	80	121880	122060	122240	122420	122600	122780	122960	123140	
	85	121890	122070	122250	122430	122610	122790	122970	123150	
	90	121900	122080	122260	122440	122620	122800	122980	123160	
	100	121910	122090	122270	122450	122630	122810	122990	123170	
	30	121920	122100	122280	122460	122640	122820	123000	123180	
	40	121930	122110	122290	122470	122650	122830	123010	123190	
	50	121940	122120	122300	122480	122660	122840	123020	123200	
10	60	121950	122130	122310	122490	122670	122850	123030	123210	
10	70	121960	122140	122320	122500	122680	122860	123040	123220	
	80	121970	122150	122330	122510	122690	122870	123050	123230	
	90	121980	122160	122340	122520	122700	122880	123060	123240	
	100	121990	122170	122350	122530	122710	122890	123070	123250	

Stirrer blade (rev.)

1	No.	No.6(rev.)	No.8(rev.)	
Outer appearance				
N	ame	Angled fan turbine	Marine blade	
Axis dia. [mm]	Blade dia. [m]	Code No.	Code No.	
	30	123260	123360	
	40	123270	123370	
	50 123280		123380	
	60	123290	123390	
8	70	123300	123400	
0	75	123310	123410	
	80	123320	123420	
	85	123330	123430	
	90	123340	123440	
	100	123350	123450	

Stirrer bar for stirrer blade X1

(Outer appearance	0
(with scr	ew)	
Product No.	0.D.×length[mm]	Material
121720	8 × 495	
121730	8×495(with screw) ※2	
121760	8 × 800	sus
121770	8 × 995	303
121780	10 × 495	
121800	10×800	
121810	10×995	

1** Does not support a stirring seal (model TBO).2** Does not support model ZZ-□□10.

Stirrer seal(vacuum stirrer) Model TBO

Code No.	Align	Outer appearance
121090	\$ 24	(A)
121100	\$ 29	
121110	\$ 34	
121120	 \$ 40	\mathbb{V}

Stirrer bar for stirrer seal

	Code No.	Material		O.D. ×length (blade size)[mm]		Oute	er appe	arance
	121520			8×450, no split 涨1		2	3	(4)
_ S	121560	SUS		8×495, with split	10°	0	n	(4) N
stirr blade	121570			8×600, with split				
ade	121490			8 × 450				
er	121500	Glass ※2		8 × 600	2			
	121510			8 × 800				
Wit	121590		T- 1	8 × 450 (80W × 17H)				
th	203210		T-11	8 × 600 (80W × 17H)			Д	Ц
st	203220	All Teflon®	T-12	8 × 700 (100W × 23H)	3			
İ	121600		T- 2	8 × 600 (120W × 23H)				
rrer	121610		T- 3	8 × 800 (150W × 25H)				
6 9	200210	All Teflon ®	T-1A	$8 \times 450 (60W \times 40H)$	4	_		
ade	200220	(anchor type)	T-2A	$8 \times 600 (80W \times 60H)$	l . A	Ø/ \	گ	

%1 In addition to a stirrer blade for stirrer seal, ordinary blades can be connected and used.

Stirrer blade for stirrer seal


Code No.	Material	Model	Standard[mm]	Outer appearance			
121450	SUS	G-1	58W×19H				
121460	303	G-2	71W×19H	1	1 2		
121430		E-1	58W×19H				
121440	Glass ※2	E-2	71W×19H				
121470	Toflon @	F-1	58W×19H	2			
121480	Teflon ®	F-2	71W×19H				

%2 When using on the model ZZ-1 10, it may break at reversal depending on the viscosity of the solvent. Stop stirring before attempting reverse operation and changing the stirring direction.

■Optional parts

P P .							
①One-touch co	upling (※3)		2)Flexible couplir	ng (※3)			
Ø							
Code No	Туре	Shaft dia.[mm]	Code No	Туре	Shaft dia.[mm]		
123480	CW-1	ϕ 6~ ϕ 10	123490	CF-1	ϕ 6~ ϕ 10		
③Flexible stirre	er (X3)		④NZ clamp				
			(Usable on 2	ZZ)	~		
				K			
(())		V			
a a					N BOOM		
					\$~		
Code No	Туре	Shaft dia.[mm]	Code No	Туре	Clamp dia.[mm]		
123570	FS-1	$\phi 1 \sim \phi 8$	203300	Model NZ	ϕ 22 (MAX)		
5Stand base	\$	~					
		> //					
	4	1 5					
		y //					
Code No	Туре	Pole size[mm]					
158710	ZS	$\phi 22 \times 780$					
203190	ZS-2	φ 22 × 1200					
		77 4 1 40					

X3 Not usable on model ZZ-1 \square 10.



①NZ arbor			②Knob bolt f	or ZZ arbor		③Plug for	NZ air purg	е
(Knob bolt	shape differ	rs from				(Same as	the accesso	rv for
that of ZZ	but can be u	sed.)	4	807		ZZ)	0	
		\supset						
				0			\mathcal{I}	
y							\checkmark	
Code No	Arbor O.D.	Qty	Code No	Standard	Qty	Code No	Tube O.D.	Qty
203270	φ16	1	256140	M6	1	203330	$\phi 6$	1
4)Stirrer bar	rscrew		⑤NZ clamp kn	ob bolt		6NZ chuck h	andle	
	ne screw inc						he accessory	for ZZ-
with stirre	r bar with s	crew)		()		1000 series	s) a	
			1	20D			$\searrow 9$	
				SAA				`
				0				/
Code No	Standard	Qty	Code No	Standard	Qty	Code No	Standard	Qty
203340	M5	1	203310	M8	1	203320		1
⑦ZZ chuck co	over		⑧Collet					
						13A		
	(PA)=B							
	The second				- 786			
	K la							
			Code No	Product	10.0100	Cumpont	d abaft	0+
	\bigcirc		Code No	Product	name		ed shaft	Qty
			256100	Collet ϕ 4		dia.	~4	1
			256100	Collet ϕ 6			~4 ~6	1
Code No	Material	Qty	256120	Collet ϕ 8			~8	1
	material	ur Ly	200120	(standard)		ψ	Ŭ	I
256150	PC resin	1	256130	Collet ϕ 10		<i>d</i> .9 <i>i</i>	~10	1
200100	1.0.100111		200100	001100910	,	φ.	. •	•